

CACNB3 Rabbit Polyclonal Antibody

Catalog No: #29659



Package Size: #29659-1 50ul #29659-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|---|
| Product Name | CACNB3 Rabbit Polyclonal Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Affinity purification |
| Applications | WB |
| Species Reactivity | Mouse,Rat |
| Immunogen Description | Recombinant fusion protein of human CACNB3 (NP_000716.2). |
| Other Names | CACNB3;CAB3;CACNLB3 |
| Accession No. | Uniprot:P54284GenelD:784 |
| Uniprot | P54284 |
| GenelD | 784 |
| Calculated MW | 50kDa/53kDa/54kDa |
| SDS-PAGE MW | 58kDa |
| Formulation | PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| Storage | Store at -20°C. Avoid freeze / thaw cycles. |

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using CACNB3 antibody.

Background

This gene encodes a regulatory beta subunit of the voltage-dependent calcium channel. Beta subunits are composed of five domains, which contribute to the regulation of surface expression and gating of calcium channels and may also play a role in the regulation of transcription factors and calcium transport.

Note: This product is for in vitro research use only