

## GPR82 Rabbit Polyclonal Antibody

Catalog No: #54232

Package Size: #54232-1 50ul #54232-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	GPR82 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Mouse,Rat
Immunogen Description	A synthetic peptide of human GPR82 (NP_543007.1).
Other Names	GPR82
Accession No.	Swiss Prot:Q96P67GeneID:27197
Uniprot	Q96P67
Calculated MW	38kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

## Application Details

WB□1:500 - 1:2000IHC□1:100 - 1:200

## Background

The protein encoded by this gene is an orphan G protein-coupled receptor of unknown function. The encoded protein is a member of a family of proteins that contain seven transmembrane domains and transduce extracellular signals through heterotrimeric G proteins.

Note: This product is for in vitro research use only