

CLCC1 Rabbit Polyclonal Antibody

Catalog No: #54448



Package Size: #54448-1 50ul #54448-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	CLCC1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human CLCC1 (NP_055942.1).
Other Names	MCLC;CLCC1
Accession No.	Swiss Prot:Q96S66GeneID:23155
Uniprot	Q96S66
Calculated MW	—
SDS-PAGE MW	62kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB□1:2000 - 1:5000IHC□1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using CLCC1 at 1:1000 dilution.

Background

—

Note: This product is for in vitro research use only