

HCAR3 Rabbit Polyclonal Antibody

Catalog No: #54531

Package Size: #54531-1 50ul #54531-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HCAR3 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human HCAR3.
Other Names	GPR109B;HCA3;HM74;PUMAG;Puma-g;HCAR3
Accession No.	Swiss Prot:P49019GenelD:8843
Uniprot	P49019
Calculated MW	—
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB□1:500 - 1:2000IHC□1:50 - 1:200

Background

—

Note: This product is for in vitro research use only