

GTF3C3 Rabbit Polyclonal Antibody

Catalog No: #54533

Package Size: #54533-1 50ul #54533-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GTF3C3 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human GTF3C3.
Other Names	TFIIIC102;TFIIICgamma;TFIIIC2-102;GTF3C3
Accession No.	Swiss Prot:Q9Y5Q9GeneID:9330
Uniprot	Q9Y5Q9
Calculated MW	—
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB□1:500 - 1:2000IHC□1:50 - 1:200IF□1:50 - 1:200

Background

—

Note: This product is for in vitro research use only