

USH1G Rabbit Polyclonal Antibody

Catalog No: #54560

Package Size: #54560-1 50ul #54560-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	USH1G Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant protein of human USH1G.
Other Names	ANKS4A;SANS;USH1G
Accession No.	Swiss Prot:Q495M9GeneID:124590
Uniprot	Q495M9
Calculated MW	—
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB□1:500 - 1:2000IF□1:50 - 1:200

Background

—

Note: This product is for in vitro research use only