

WBP2NL Rabbit Polyclonal Antibody

Catalog No: #54567

Package Size: #54567-1 50ul #54567-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	WBP2NL Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant protein of human WBP2NL.
Other Names	GRAMD7;PAWP;WBP2NL
Accession No.	Swiss Prot:Q6ICG8GeneID:164684
Uniprot	Q6ICG8
Calculated MW	—
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB□1:500 - 1:2000IF□1:50 - 1:200

Background

—

Note: This product is for in vitro research use only