

# DGCR8 Rabbit Polyclonal Antibody

Catalog No: #54579



Package Size: #54579-1 50ul #54579-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                         |
|-----------------------|---------------------------------------------------------|
| Product Name          | DGCR8 Rabbit Polyclonal Antibody                        |
| Host Species          | Rabbit                                                  |
| Clonality             | Polyclonal                                              |
| Isotype               | IgG                                                     |
| Purification          | Affinity purification                                   |
| Applications          | WB,IHC                                                  |
| Species Reactivity    | Human,Mouse,Rat                                         |
| Immunogen Description | Recombinant protein of human DGCR8.                     |
| Other Names           | —                                                       |
| Accession No.         | Swiss Prot:Q8WYQ5GeneID:54487                           |
| Uniprot               | Q8WYQ5                                                  |
| Calculated MW         | —                                                       |
| SDS-PAGE MW           | 110kDa                                                  |
| Formulation           | Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| Storage               | Store at -20°C. Avoid freeze / thaw cycles.             |

## Application Details

WB 1:2000 - 1:5000IHC 1:50 - 1:200

## Images



Western blot analysis of extracts of various cell lines, using DGCR8 at 1:1000 dilution.

## Background

—

---

Note: This product is for in vitro research use only