

SIGLEC15 Rabbit Polyclonal Antibody

Catalog No: #54584

Package Size: #54584-1 50ul #54584-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SIGLEC15 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Rat
Immunogen Description	Recombinant protein of human Siglec-15.
Other Names	CD33L3, HsT1361, SIGLEC-15
Accession No.	Uniprot:Q6ZMC9GeneID:284266
Uniprot	Q6ZMC9
GeneID	284266
Calculated MW	—
SDS-PAGE MW	36KDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200IF 1:50 - 1:200


Images




Immunohistochemistry of paraffin-embedded human spleen using Siglec-15 Rabbit pAb.




Immunohistochemistry of paraffin-embedded human esophageal cancer using Siglec-15 Rabbit pAb.



Immunofluorescence analysis of Rat spleen cells using Siglec-15 antibody.



Immunofluorescence analysis of Human spleen cells using Siglec-15 antibody.



Western blot analysis of extracts of various cell lines, using SIGLEC15 antibody.

Background

SIGLEC15 (Sialic Acid Binding Ig Like Lectin 15) is a Protein Coding gene. Diseases associated with SIGLEC15 include Osteoporosis, Juvenile and Osteoporosis. Among its related pathways are Innate Immune System and RET signaling. An important paralog of this gene is SIGLEC1.

Note: This product is for in vitro research use only