

NUB1 Rabbit Polyclonal Antibody

Catalog No: #54585



Package Size: #54585-1 50ul #54585-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NUB1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human IFT140.
Other Names	—
Accession No.	Swiss Prot:Q9Y5A7GeneID:51667
Uniprot	Q9Y5A7
Calculated MW	—
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB□1:500 - 1:2000IF□1:50 - 1:200

Background

—

Note: This product is for in vitro research use only