

MAPK8/MAPK9/MAPK10 Rabbit Polyclonal Antibody

Catalog No: #54589



Package Size: #54589-1 50ul #54589-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	MAPK8/MAPK9/MAPK10 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human MAPK8/MAPK9/MAPK10.
Other Names	—
Accession No.	Swiss Prot:P45983/P45984/P53779GenelD:5599/5601/5602
Uniprot	P45983/P45984/P53779
Calculated MW	—
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB□1:2000 - 1:5000IHC□1:50 - 1:200

Background

—

Note: This product is for in vitro research use only