

CUL2 Rabbit Polyclonal Antibody

Catalog No: #54727

Package Size: #54727-1 50ul #54727-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CUL2 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human CUL2 (NP_001185708.1).
Other Names	CUL2;cullin-2
Accession No.	Swiss Prot:Q13617GeneID:8453
Uniprot	Q13617
Calculated MW	86kDa/89kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB□1:500 - 1:2000IHC□1:50 - 1:200IF□1:50 - 1:200

Images



Immunohistochemistry of paraffin-embedded mouse heart using CUL2 at dilution of 1:200 (40x lens).



Immunofluorescence analysis of MCF-7 cells using CUL2 .
Blue: DAPI for nuclear staining.

Background

-

Note: This product is for in vitro research use only