

KIFC1 Rabbit Polyclonal Antibody

Catalog No: #54754

Package Size: #54754-1 50ul #54754-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	KIFC1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human KIFC1 (NP_002254.2).
Other Names	KIFC1;HSET;KNSL2
Accession No.	Swiss Prot:Q9BW19GeneID:3833
Uniprot	Q9BW19
Calculated MW	73kDa
SDS-PAGE MW	74kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:100IF 1:50 - 1:100


Images




Immunohistochemistry of paraffin-embedded rat testis using KIFC1 at dilution of 1:100 (40x lens).




Immunohistochemistry of paraffin-embedded human colon carcinoma using KIFC1 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse testis using KIFC1 at dilution of 1:100 (40x lens).



Western blot analysis of extracts of various cell lines, using KIFC1 at 1:3000 dilution.



Immunofluorescence analysis of U-2 OS cells using KIFC1 Polyclonal at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.

Background

Note: This product is for in vitro research use only