

CDK14 Rabbit Polyclonal Antibody

Catalog No: #54758



Package Size: #54758-1 50ul #54758-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	CDK14 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human CDK14
Other Names	CDK14;PFTAIR1;PFTK1
Accession No.	Swiss Prot:O94921GeneID:5218
Uniprot	O94921
Calculated MW	47kDa/50kDa/53kDa
SDS-PAGE MW	45kDa
Formulation	Buffer: PBS with 0.02% sodium azide, pH7.3.
Storage	Store at 4°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:1000IHC 1:50 - 1:100FC 1:20 - 1:50

Images



Western blot analysis of extracts of human brain, using CDK14 .

Background

-

Note: This product is for in vitro research use only