

UQCRC1 Rabbit Polyclonal Antibody

Catalog No: #54759

Package Size: #54759-1 50ul #54759-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	UQCRC1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human UQCRC1 (NP_003356.2).
Other Names	UQCRC1;D3S3191;QCR1;UQCR1
Accession No.	Swiss Prot:P31930GenelD:7384
Uniprot	P31930
Calculated MW	52kDa
SDS-PAGE MW	53kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200


Images




Immunohistochemistry of paraffin-embedded rat brain using UQCRC1 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human prostate using UQCRC1 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human stomach using UQCRC1 at dilution of 1:100 (40x lens).



Western blot analysis of extracts of various cell lines, using UQCRC1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only