

ARHGAP5 Rabbit Polyclonal Antibody

Catalog No: #54774

Package Size: #54774-1 50ul #54774-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ARHGAP5 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human RhoGAP (NP_001164.2).
Other Names	ARHGAP5;GFI2;RhoGAP5;p190-B;p190BRhoGAP
Accession No.	Uniprot:Q13017GeneID:394
Uniprot	Q13017
GeneID	394
Calculated MW	26kDa/27kDa/172kDa
SDS-PAGE MW	185kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB□1:500 - 1:2000IF□1:50 - 1:200

Images



Immunofluorescence analysis of HeLa cells using RhoGAP antibody.



Western blot analysis of extracts of various cell lines, using RhoGAP antibody.

Background

Rho GTPase activating protein 5 negatively regulates RHO GTPases, a family which may mediate cytoskeleton changes by stimulating the hydrolysis of bound GTP. Two transcript variants encoding different isoforms have been found for this gene.

Note: This product is for in vitro research use only