

HMGCS1 Rabbit Polyclonal Antibody

Catalog No: #54807



Package Size: #54807-1 50ul #54807-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	HMGCS1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human HMGCS1 (NP_002121.4).
Other Names	HMGCS1;HMGCS
Accession No.	Swiss Prot:Q01581GeneID:3157
Uniprot	Q01581
Calculated MW	57kDa
SDS-PAGE MW	57kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:1000 - 1:2000IF 1:50 - 1:200


Images



Western blot analysis of extracts of various cell lines, using HMGCS1 at 1:7000 dilution.



Immunofluorescence analysis of C6 cells using HMGCS1 Polyclonal at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.



Immunofluorescence analysis of HeLa cells using HMGCS1 Polyclonal at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.

Background

—

Note: This product is for in vitro research use only