

ABCF2 Rabbit Polyclonal Antibody

Catalog No: #54851



Package Size: #54851-1 50ul #54851-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ABCF2 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human ABCF2 (NP_005683.2).
Other Names	ABCF2;ABC28;EST133090;HUSSY-18;HUSSY18
Accession No.	Uniprot:Q9UG63GeneID:10061
Uniprot	Q9UG63
GeneID	10061
Calculated MW	71kDa/72kDa
SDS-PAGE MW	72kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using ABCF2 antibody.



Immunofluorescence analysis of U2OS cells using ABCF2 antibody.

Background

This gene encodes a member of the ATP-binding cassette (ABC) transporter superfamily. ATP-binding cassette proteins transport various molecules across extra- and intracellular membranes. Alterations in this gene may be involved in cancer progression. Alternative splicing results in multiple transcript variants. Related pseudogenes have been identified on chromosomes 3 and 7.

Note: This product is for in vitro research use only