

BCAS2 Rabbit Polyclonal Antibody

Catalog No: #54853



Package Size: #54853-1 50ul #54853-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	BCAS2 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human BCAS2 (NP_005863.1).
Other Names	BCAS2;DAM1;SPF27;Snt309
Accession No.	Swiss Prot:O75934GeneID:10286
Uniprot	O75934
Calculated MW	26kDa
SDS-PAGE MW	25kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200IF 1:50 - 1:200


Images




Western blot analysis of extracts of various cell lines, using BCAS2 at 1:1000 dilution.




Immunohistochemistry of paraffin-embedded human lung cancer using BCAS2 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human esophagus using BCAS2 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse brain using BCAS2 at dilution of 1:100 (40x lens).



Confocal immunofluorescence analysis of U2OS cells using BCAS2 Polyclonal at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

—

Note: This product is for in vitro research use only