

CSTF2T Rabbit Polyclonal Antibody

Catalog No: #54867



Package Size: #54867-1 50ul #54867-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	CSTF2T Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human CSTF2T (NP_056050.1).
Other Names	CSTF2T;CstF-64T
Accession No.	Swiss Prot:Q9H0L4GeneID:23283
Uniprot	Q9H0L4
Calculated MW	64kDa
SDS-PAGE MW	70kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200IF 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using CSTF2T at 1:1000 dilution.

Background

—

Note: This product is for in vitro research use only