

CBX6 Rabbit Polyclonal Antibody

Catalog No: #54913



Package Size: #54913-1 50ul #54913-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	CBX6 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human CBX6 (NP_055107.3).
Other Names	CBX6
Accession No.	Swiss Prot:O95503GeneID:23466
Uniprot	O95503
Calculated MW	43kDa
SDS-PAGE MW	50kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200


Images




Western blot analysis of extracts of HeLa cells, using CBX6 at 1:1000 dilution.



Immunohistochemistry of paraffin-embedded rat liver using CBX6 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded rat brain using CBX6 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse testis using CBX6 at dilution of 1:100 (40x lens).

Background

Note: This product is for in vitro research use only