

ATG9A Rabbit Polyclonal Antibody

Catalog No: #54918

Package Size: #54918-1 50ul #54918-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ATG9A Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human ATG9A (NP_076990.4).
Other Names	ATG9A;APG9L1;MGD3208;mATG9
Accession No.	Swiss Prot:Q7Z3C6GeneID:79065
Uniprot	Q7Z3C6
Calculated MW	60kDa/87kDa/94kDa
SDS-PAGE MW	100kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:1000 - 1:2000IF 1:50 - 1:200

Images



Immunofluorescence analysis of U2OS cells using ATG9A at dilution of 1:100. Blue: DAPI for nuclear staining.



Western blot analysis of extracts of 293T cells, using ATG9A .

Background

—

Note: This product is for in vitro research use only