

SLC8A1 Rabbit Polyclonal Antibody

Catalog No: #54919



Package Size: #54919-1 50ul #54919-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	SLC8A1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human SLC8A1 (NP_066920.1).
Other Names	SLC8A1;NCX1
Accession No.	Swiss Prot:P32418GenelD:6546
Uniprot	P32418
Calculated MW	104kDa/107kDa/108kDa
SDS-PAGE MW	109kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:2000 - 1:5000IHC 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using SLC8A1 at 1:1000 dilution.

Background

-

Note: This product is for in vitro research use only