

IRAK2 Rabbit Polyclonal Antibody

Catalog No: #55005



Package Size: #55005-1 50ul #55005-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	IRAK2 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human IRAK2 (NP_001561.3).
Other Names	IRAK2;IRAK-2
Accession No.	Uniprot:O43187GeneID:3656
Uniprot	O43187
GeneID	3656
Calculated MW	69kDa
SDS-PAGE MW	75kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using IRAK2 antibody.

Background

IRAK2 encodes the interleukin-1 receptor-associated kinase 2, one of two putative serine/threonine kinases that become associated with the interleukin-1 receptor (IL1R) upon stimulation. IRAK2 is reported to participate in the IL1-induced upregulation of NF-kappaB.

Note: This product is for in vitro research use only