

TRMT1 Rabbit Polyclonal Antibody

Catalog No: #55126



Package Size: #55126-1 50ul #55126-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	TRMT1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human TRMT1 (NP_060192.1).
Other Names	TRMT1;TRM1
Accession No.	Uniprot:Q9NXH9GeneID:55621
Uniprot	Q9NXH9
GeneID	55621
Calculated MW	69kDa/72kDa
SDS-PAGE MW	90kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using TRMT1 antibody.

Background

Dimethylates a single guanine residue at position 26 of most tRNAs using S-adenosyl-L-methionine as donor of the methyl groups.

Note: This product is for in vitro research use only