

ZFAND3 Rabbit Polyclonal Antibody

Catalog No: #55228



Package Size: #55228-1 50ul #55228-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ZFAND3 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ZFAND3 (NP_068762.1).
Other Names	ZFAND3;TEX27
Accession No.	Swiss Prot:Q9H8U3GeneID:60685
Uniprot	Q9H8U3
Calculated MW	25kDa
SDS-PAGE MW	27kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200


Images



Western blot analysis of extracts of various cell lines, using ZFAND3 at 1:1000 dilution.



Immunohistochemistry of paraffin-embedded rat kidney using ZFAND3 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse kidney using ZFAND3 at dilution of 1:100 (40x lens).

Background

Note: This product is for in vitro research use only