

STRN Rabbit Polyclonal Antibody

Catalog No: #55305



Package Size: #55305-1 50ul #55305-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	STRN Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human STRN (NP_003153.2).
Other Names	STRN;PPP2R6A;SG2NA;STRN1;striatin
Accession No.	Uniprot:O43815GeneID:6801
Uniprot	O43815
GeneID	6801
Calculated MW	80kDa/86kDa
SDS-PAGE MW	120kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB1:500 - 1:2000IHC1:50 - 1:200


Images




Western blot analysis of extracts of various cell lines, using STRN antibody.



Immunohistochemistry of paraffin-embedded human liver damage using STRN antibody.



Immunohistochemistry of paraffin-embedded rat heart using STRN antibody.



Immunohistochemistry of paraffin-embedded mouse brain using STRN antibody.

Background

Calmodulin-binding protein which may function as scaffolding or signaling protein and may play a role in dendritic Ca^{2+} signaling.

Note: This product is for in vitro research use only