

CKAP4 Rabbit Polyclonal Antibody

Catalog No: #55326



Package Size: #55326-1 50ul #55326-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	CKAP4 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human CKAP4
Other Names	CKAP4;CLIMP-63;ERGIC-63;p63
Accession No.	Swiss Prot:Q07065GeneID:10970
Uniprot	Q07065
Calculated MW	66kDa
SDS-PAGE MW	66kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:100IF 1:50 - 1:100


Images




Western blot analysis of extracts of various cell lines, using CKAP4 at 1:500 dilution.




Immunohistochemistry of paraffin-embedded mouse brain using CKAP4 (40x lens).



Immunohistochemistry of paraffin-embedded rat brain using CKAP4 (40x lens).



Immunohistochemistry of paraffin-embedded mouse lung using CKAP4 (40x lens).



Confocal immunofluorescence analysis of U2OS cells using CKAP4 Polyclonal at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

—

Note: This product is for in vitro research use only