

CEACAM7 Rabbit Polyclonal Antibody

Catalog No: #55427



Package Size: #55427-1 50ul #55427-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	CEACAM7 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human CEACAM7 (NP_008821.2).
Other Names	CEACAM7;CGM2
Accession No.	Uniprot:Q14002GeneID:1087
Uniprot	Q14002
GeneID	1087
Calculated MW	19kDa/29kDa
SDS-PAGE MW	29kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using CEACAM7 antibody.

Background

This gene encodes a cell surface glycoprotein and member of the carcinoembryonic antigen (CEA) family of proteins. Expression of this gene may be downregulated in colon and rectal cancer. Additionally, lower expression levels of this gene may be predictive of rectal cancer recurrence. This gene is present in a CEA family gene cluster on chromosome 19. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only