

UBE2U Rabbit Polyclonal Antibody

Catalog No: #55494



Package Size: #55494-1 50ul #55494-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	UBE2U Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human UBE2U (NP_689702.1).
Other Names	UBE2U
Accession No.	Swiss Prot:Q5VVX9GeneID:148581
Uniprot	Q5VVX9
Calculated MW	26kDa/37kDa
SDS-PAGE MW	42kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:100


Images



Western blot analysis of extracts of various cell lines, using UBE2U at 1:1000 dilution.



Immunohistochemistry of paraffin-embedded human liver using UBE2U at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse brain using UBE2U at dilution of 1:100 (40x lens).

Background

Note: This product is for in vitro research use only