

MARK1 Rabbit Polyclonal Antibody

Catalog No: #55514



Package Size: #55514-1 50ul #55514-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	MARK1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human MARK1 (NP_061120.3).
Other Names	MARK1;MARK;Par-1c;Par1c
Accession No.	Swiss Prot:Q9P0L2GeneID:4139
Uniprot	Q9P0L2
Calculated MW	72kDa/84kDa/89kDa
SDS-PAGE MW	95kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using MARK1 at 1:1000 dilution.



Immunohistochemistry of paraffin-embedded mouse lung using MARK1 at dilution of 1:100 (40x lens).

Background

-

Note: This product is for in vitro research use only