

MCTS1 Rabbit Polyclonal Antibody

Catalog No: #29287



Package Size: #29287-1 50ul #29287-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	MCTS1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human MCTS1 (NP_054779.1).
Other Names	MCTS1;MCT-1;MCT1
Accession No.	Swiss Prot:Q9ULC4GeneID:28985
Calculated MW	19kDa/20kDa
SDS-PAGE MW	20kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details


WB 1:500 - 1:2000

IF 1:50 - 1:200


Images




Western blot analysis of extracts of various cell lines, using MCTS1 at 1:3000 dilution.



Immunofluorescence analysis of C6 cells using MCTS1 at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.



Immunofluorescence analysis of NIH/3T3 cells using MCTS1 at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.



Immunofluorescence analysis of U-2OS cells using MCTS1 at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.

Background

—

Note: This product is for in vitro research use only