

PMAIP1 Rabbit Polyclonal Antibody

Catalog No: #55610



Package Size: #55610-1 50ul #55610-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PMAIP1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human PMAIP1 (NP_066950.1).
Other Names	PMAIP1;APR;NOXA
Accession No.	Swiss Prot:Q13794GeneID:5366
Uniprot	Q13794
Calculated MW	6kDa/14kDa
SDS-PAGE MW	14kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200


Images




Western blot analysis of extracts of various cell lines, using PMAIP1 at 1:1000 dilution.




Immunofluorescence analysis of Jurkat cells using PMAIP1 at dilution of 1:100. Blue: DAPI for nuclear staining.




Immunofluorescence analysis of rat breast cells using PMAIP1 at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of mouse breast cells using PMAIP1 at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of mouse breast cells using PMAIP1 at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of mouse breast cells using PMAIP1 at dilution of 1:100. Blue: DAPI for nuclear staining.

—

Note: This product is for in vitro research use only