

Rabbit Myoglobin (MYO) ELISA Kit

Catalog No: #EK20415



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Rabbit Myoglobin (MYO) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rabbit
Storage	12 months

Application Details

Detect Range:78.1-5000pg/mL
Sensitivity:28pg/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only