

## Rat Interleukin 10 (IL10) ELISA Kit

Catalog No: #EK20209

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                     |
|--------------------|-------------------------------------|
| Product Name       | Rat Interleukin 10 (IL10) ELISA Kit |
| Brief Description  | ELISA Kit                           |
| Applications       | ELISA                               |
| Species Reactivity | Rat                                 |
| Storage            | 12 months                           |

## Application Details

Detect Range:78.1-5000pg/mL

Sensitivity:27pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &amp;mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only