

Rat Tumor Necrosis Factor Alpha (TNFa) ELISA Kit

Catalog No: #EK20211

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Tumor Necrosis Factor Alpha (TNFa) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rat
Storage	12 months

Application Details

Detect Range:15.6-1000pg/ml

Sensitivity:7.1pg/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only