

TAZ (phospho-Ser89) rabbit pAb

Catalog No: #13489

Package Size: #13489-1 50ul #13489-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TAZ (phospho-Ser89) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	This antibody detects endogenous levels of Human Mouse Rat TAZ (phospho-Ser89)
Immunogen Description	Synthesized phosho peptide around human TAZ (Ser89)
Other Names	WW domain-containing transcription regulator protein 1 (Transcriptional coactivator with PDZ-binding motif)
Accession No.	Swiss Prot:Q9GZV5GeneID:25937
Uniprot	Q9GZV5
GeneID	25937
SDS-PAGE MW	44
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only