

Scribble (phospho-Ser1220) rabbit pAb

Catalog No: #13536



Package Size: #13536-1 50ul #13536-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Scribble (phospho-Ser1220) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human, Mouse
Specificity	This antibody detects endogenous levels of Human Mouse Scribble (phospho-Ser1220)
Immunogen Description	Synthesized phospho peptide around human Scribble (Ser1220)
Other Names	Protein scribble homolog (Scribble) (hScrib) (Protein LAP4)
Accession No.	Swiss Prot:Q14160 GeneID:23513
Uniprot	Q14160
GeneID	23513
SDS-PAGE MW	179
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Background

scribbled planar cell polarity protein (SCRIB) Homo sapiens This gene encodes a protein that was identified as being similar to the Drosophila scribble protein. The mammalian protein is involved in tumor suppression pathways. As a scaffold protein involved in cell polarization processes, this protein binds to many other proteins. The encoded protein binds to papillomavirus E6 protein via its PDZ domain and the C-terminus of E6. Two alternatively spliced transcript variants that encode different protein isoforms have been found for this gene. [provided by RefSeq, Nov 2011],

Note: This product is for in vitro research use only