

RIP (phospho-Ser166) rabbit pAb

Catalog No: #13557



Package Size: #13557-1 50ul #13557-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RIP (phospho-Ser166) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	This antibody detects endogenous levels of Human RIP (phospho-Ser166)
Immunogen Description	Synthesized phosho peptide around human RIP (Ser166)
Other Names	Receptor-interacting serine/threonine-protein kinase 1 (EC 2.7.11.1) (Cell death protein RIP) (Receptor-interacting protein 1) (RIP-1) (Serine/threonine-protein kinase RIP)
Accession No.	Swiss Prot:Q13546GeneID:8737
Uniprot	Q13546
GeneID	8737
SDS-PAGE MW	73
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only