

## MERIT40 (phospho-Ser29) rabbit pAb

Catalog No: #13710



Package Size: #13710-1 50ul #13710-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	MERIT40 (phospho-Ser29) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of Human MERIT40 (phospho-Ser29)
Immunogen Description	Synthesized phospho peptide around human MERIT40 (Ser29)
Other Names	BRISC and BRCA1-A complex member 1 (Mediator of RAP80 interactions and targeting subunit of 40 kDa) (New component of the BRCA1-A complex)
Accession No.	Swiss Prot:Q9NWW8GenelD:29086
Uniprot	Q9NWW8
GenelD	29086
SDS-PAGE MW	36
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

## Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only