

Keratin 17 (phospho-Ser44) rabbit pAb

Catalog No: #13766

Package Size: #13766-1 50ul #13766-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Keratin 17 (phospho-Ser44) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of Human Keratin 17 (phospho-Ser44)
Immunogen Description	Synthesized phosho peptide around human Keratin 17 (Ser44)
Other Names	Keratin, type I cytoskeletal 17 (39.1) (Cytokeratin-17) (CK-17) (Keratin-17) (K17)
Accession No.	Swiss Prot:Q04695GeneID:3872
Uniprot	Q04695
GeneID	3872
SDS-PAGE MW	48
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Background

keratin 17(KRT17) Homo sapiens This gene encodes the type I intermediate filament chain keratin 17, expressed in nail bed, hair follicle, sebaceous glands, and other epidermal appendages. Mutations in this gene lead to Jackson-Lawler type pachyonychia congenita and steatocystoma multiplex. [provided by RefSeq, Aug 2008].

Note: This product is for in vitro research use only