

JunB (phospho-Thr102/Thr104) rabbit pAb

Catalog No: #13768



Package Size: #13768-1 50ul #13768-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	JunB (phospho-Thr102/Thr104) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of Human JunB (phospho-Thr102 or Thr104)
Immunogen Description	Synthesized phospho peptide around human JunB (Thr102 and Thr104)
Other Names	Transcription factor jun-B
Accession No.	Swiss Prot:P17275GeneID:3726
Uniprot	P17275
GeneID	3726
SDS-PAGE MW	36
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only