

## FoxO4 (phospho-Ser193) rabbit pAb

Catalog No: #13847

Package Size: #13847-1 50ul #13847-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	FoxO4 (phospho-Ser193) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of Human FoxO4 (phospho-Ser193)
Immunogen Description	Synthesized phospho peptide around human FoxO4 (Ser193)
Other Names	Forkhead box protein O4 (Fork head domain transcription factor AFX1)
Accession No.	Swiss Prot:P98177GeneID:4303
Uniprot	P98177
GeneID	4303
SDS-PAGE MW	55
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

## Application Details

WB 1:1000-2000

## Background

forkhead box O4(FOXO4) Homo sapiens This gene encodes a member of the O class of winged helix/forkhead transcription factor family. Proteins encoded by this class are regulated by factors involved in growth and differentiation indicating they play a role in these processes. A translocation involving this gene on chromosome X and the homolog of the Drosophila trithorax gene, encoding a DNA binding protein, located on chromosome 11 is associated with leukemia. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2010],

Note: This product is for in vitro research use only