

BCL9L (phospho-Ser915) rabbit pAb

Catalog No: #14050

Package Size: #14050-1 50ul #14050-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	BCL9L (phospho-Ser915) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of Human BCL9L (phospho-Ser915)
Immunogen Description	Synthesized phosho peptide around human BCL9L (Ser915)
Other Names	B-cell CLL/lymphoma 9-like protein (B-cell lymphoma 9-like protein) (BCL9-like protein) (Protein BCL9-2)
Accession No.	Swiss Prot:Q86UU0GeneID:283149
Uniprot	Q86UU0
GeneID	283149
SDS-PAGE MW	160
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only