

## SARS-CoV-2 (2019-nCoV) S1 protein NTD, hFc Tag

Catalog No: #AP89616

Package Size: #AP89616-1 10ug #AP89616-2 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	SARS-CoV-2 (2019-nCoV) S1 protein NTD, hFc Tag
Host Species	HEK293
Purification	The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.
Species Reactivity	Human
Immunogen Description	S1 protein NTD(Ser13-Leu303)+hFc(Glu99-Ala330)
Other Names	S1 protein NTD,Spike protein S1 NTD,BetaCoV S1-NTD,COVID-19
Calculated MW	59.0 kDa
Tag Info	C-Human Fc tag
Formulation	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 96 % trehalose is added as protectants before lyophilization.
Storage	Store at -80°C for 12 months (Avoid repeated freezing and thawing)

## Images





Figure 1. SARS-CoV-2 (2019-nCoV) S1 protein NTD, hFc Tag on SDS-PAGE under reducing condition.



## Product Description

Expression Region:550Research Topic:SARS-CoV-2 (Severe Acute Respiratory Syndrome Coronavirus 2) also known as Covid19 (2019 Novel Coronavirus) is a virus that causes illnesses ranging from the common cold to severe diseases. The spike protein is a type I transmembrane protein

containing two subunits, S1 and S2. S1 mainly contains a receptor binding domain (RBD), which accounts for recognizing the cell surface receptor, ACE2. S2 contains basic elements needed for the membrane fusion. Recent publications indicate that S1-RBD domain can induce virus neutralizing-antibody and T cell response.

---

Note: This product is for in vitro research use only