

## BCLAF1 (Phospho-Ser512) Antibody

Catalog No: #SAB633

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	BCLAF1 (Phospho-Ser512) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified from rabbit serum by affinity purification via sequential chromatography on phospho-peptide and non-phospho-peptide affinity columns.
Applications	WB
Species Reactivity	HM
Specificity	BCLAF1(Phospho-Ser512) Antibody detects endogenous levels of BCLAF1 only when phosphorylated at Ser512.
Immunogen Description	A synthesized peptide derived from human BCLAF1 around the phosphorylation site of Ser512.
Other Names	Bcl-2-associated transcription factor 1; Btf; BCLAF1 and THRAP3 family member 1
Uniprot	Q9NYF8
GeneID	9774
SDS-PAGE MW	106kDa
Concentration	1 mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C/1 year

## Application Details

Western Blot: 1/500 - 1/2000

## Background

Death-promoting transcriptional repressor. May be involved in cyclin-D1/CCND1 mRNA stability through the SNARP complex which associates with both the 3'end of the CCND1 gene and its mRNA.

Note: This product is for in vitro research use only