

Rat Carbohydrate Antigen 125 (CA125) ELISA Kit

Catalog No: #EK20147



Package Size: #EK20147-1 48T #EK20147-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Carbohydrate Antigen 125 (CA125) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rt

Application Details

Detect Range:1.56-100U/ml

Sensitivity:0.55U/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only