

Recombinant 2019-nCoV S Protein RBD (K417N, E484K, N501Y, Mammalian, C-His)


Catalog No: #AG0071

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant 2019-nCoV S Protein RBD (K417N, E484K, N501Y, Mammalian, C-His)
Purification	>90% as determined by SDS-PAGE
Other Names	RBD Protein, Spike RBD Protein, South African (ZA) variant (B.1.351), Beta
Calculated MW	35kDa
Formulation	Supplied as solution form in PBS, pH7.5 or lyophilized from PBS, pH7.5
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week .Store at -20 to -80 °C for twelve months from the date of receipt.

Images



Product Description

In general, proteins are shipped out with blue ice unless customers require otherwise.

Note: This product is for in vitro research use only