

MOAP1 Polyclonal Antibody

Catalog No: #29571

Package Size: #29571-1 50ul #29571-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MOAP1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Synthesized peptide derived from part region of human protein
Calculated MW	38kDa
Concentration	1mg/ml
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000

Images



Western blot analysis of lysates from A431 cells, primary antibody was diluted at 1:1000, 4° over night

Background

modulator of apoptosis 1(MOAP1) Homo sapiens The protein encoded by this gene was identified by its interaction with apoptosis regulator BAX protein. This protein contains a Bcl-2 homology 3 (BH3)-like motif, which is required for the association with BAX. When overexpressed, this gene has been shown to mediate caspase-dependent apoptosis.

Note: This product is for in vitro research use only