

Recombinant Human C-X-C motif chemokine 10 protein(CXCL10)

Catalog No: #AP73359

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #AP73359-1 5ug #AP73359-2 25ug #AP73359-3 100ug #AP73359-4 250ug

#AP73359-5 500ug #AP73359-6 1mg

Description

Product Name	Recombinant Human C-X-C motif chemokine 10 protein(CXCL10)
Brief Description	Recombinant Protein
Host Species	E.Coli
Purification	>97% as determined by SDS-PAGE and HPLC.
Immunogen Description	Expression Region:22-98aaSequence Info:Full Length of Mature Protein
Other Names	10 kDa interferon gamma-induced protein, IP-10, Small-inducible cytokine B10, .]
Accession No.	P02778
Uniprot	P02778
GeneID	3627;
Calculated MW	8.6 kDa
Tag Info	NO-tagged
Target Sequence	VPLSRTVRCT CISISNPVN PRSLEKLEII PASQFCPRVE IIATMKKKGE KRCLNPESKA IKNLLKAVSK EMSKRSP
Formulation	0.2 ?m filtered PBS, pH 7.4 ,lyophilized
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Chemotactic for monocytes and T-lymphocytes. Binds to CXCR3.

References

PubMed ID: 3925348; 15489334; 8423327; 10233762; 11559369; 18645041; 12173928Research Topic:Immunology

Note: This product is for in vitro research use only